

*TLink*TM DATA MANAGEMENT SYSTEM

An automated data acquisition system that allows primary focus on the patient



TLink™ Data Management System is a data acquisition system so automated it allows the perfusionist's primary focus to be on the patient, not on the screen. It incorporates all the functionality needed to implement a paperless system for perfusion records. Loaded with features, it can be integrated into a perfusion practice at a comfortable pace.

Customizable


TLink System's settings and features can be tailored to meet the specific requirements of different cases and different users, then saved and recalled as needed.

- Select the features you want to use.
- Select the data to collect, how often to collect it, and how to view it.
- Change the configuration as your practice changes.
- Manage the data to suit your requirements.
- Data entered from devices without user intervention.
- Customized screen layouts provide unlimited pump record designs.
- Both touch screen and laptop configurations available.
- User-selected security levels.
- Easily establish quality assurance attributes for continuous quality improvement (CQI) management.



Easy to use

- TLink System is an easy solution to paperless records.
- The intuitive interface, drop-down menus and function keys are easy to use for fast selection of events. One or two clicks and it's done.
- Easy barcode scanning.
- Data is acquired quickly and accurately without manual intervention, so the patient remains the perfusionist's primary focus.



Connectivity – Automatic data acquisition

- Advanced connectivity to devices in the operating room environment. Contact Terumo Cardiovascular Systems for the most up-to-date list of devices that interface with the TLink System.



Built-in Query Manager

TLink System's query management is integrated into the system.

- Easy to learn and use.
- All data can be queried within the system.
- Enhanced open-ended query capabilities can optimize your analysis post case.
- A robust database stores data efficiently and provides easy manipulation of data whenever necessary.
- Search patient records for information relevant to studies.

Flexible reporting capability

TLink System uses a variety of reports, including:

- Perfusion Case Report
- Clinical Activity
- Perfusion Checklist
- Quality Assurance
- Audit Summary Report
- Audit Detail Report

Real system support


Terumo provides its full commitment and a real support system.

Managing case data


Your data collection and storage is organized and accurate with the TLink System.

- Easily transfer cases to a hospital information system (HIS).
- Centralize all case data.
- Quick data transfer via wireless, CD, LAN or USB device.
- An easy solution to paperless records.
- Query capabilities within the system.


Bracket options provide greater flexibility



Notebook small tray bracket, 815689



Touch screen small tray bracket, 814806 and extension pole bracket, 814858



Touch screen mounting bracket, 814857

TLink Data Management System

Ordering Information

Description	Catalog No.
TLink™ Data Management System	
TLink™ Data Management System	814850
TLink™ Data Management software license	814871
Accessories	
TLink™ System touch screen computer	814851
Notebook computer	814852
Barcode scanner	814869
Touch screen mounting bracket	814857
Touch screen small tray bracket	814806
Touch screen large tray bracket	814855
Notebook small tray bracket	815689
Notebook large tray bracket	814856
Extension pole mount bracket	814858
Serial converter hub	814854
Single port serial converter	814808
Touch screen stylus	814870

Minimum hardware requirements

- Processor: Pentium 4 Class running at 1.8GHz
- Memory: 1 GB
- Display: Resolution 1280 horizontal
- Hard Disk Drive: 10 GB of available hard disk space
- CD Drive: 24X CDRW
- USB: 2 ports
- Network: PCMCIA slot available or wireless network ready

TLink™ Data Management System



Manufactured by Terumo Cardiovascular Systems, Ann Arbor, Michigan, USA.

For more information on this or other Terumo cardiovascular products, contact:

TERUMO CARDIOVASCULAR SYSTEMS CORPORATION
6200 Jackson Road
Ann Arbor, Michigan 48103-9300
USA
734 663 4145 phone
734 663 7981 fax
800 521 2818 toll free
www.terumo-cvs.com

TERUMO CORPORATION
44-1, 2-chome
Hatagaya, Shibuya-ku
Tokyo 151-0072
Japan
81 3 3374 8111 phone
81 3 3374 8196 fax
www.terumo.com

TERUMO EUROPE N.V.
Authorized EC Representative
Interleuvenlaan 40 B-3001
Leuven
Belgium
32 16 38 12 11 phone
32 16 40 02 49 fax
www.terumo-europe.com

TERUMO EUROPE N.V.
Cardiovascular Division
Hauptstrasse 87
D-65760 Eschborn
Germany
49 6196 8023 500 phone
49 6196 8023 555 fax

TERUMO LATIN AMERICA CORPORATION
8750 NW 36th Street, Suite 600
Miami, Florida 33178
USA
305 477 4822 phone
305 477 4872 fax
800 283 7866 toll free